



## Woven Hinged Knee Support

**Ref: 116**

Manufactured of flexible and sweat-preventing woven fabric. Inner silicon pad support patella perimeter. The underwire on both sides support knee joint, reduce pain by relieving the load on the joint. Velcro fasteners on the upper and lower calves increase stabilization.

### Indications:

- 1-** Ligament injury
- 2-** Meniscus rapture
- 3-** Osteoarthritis, bursitis, tendinitis
- 4-** Post-op or trauma for protection and rehabilitation purposes

### Size;

- S** : 31 - 34 cm
- M** : 34 - 37 cm
- L** : 37 - 40 cm
- XL** : 40 - 44 cm
- XXL** : 44 - 48 cm

**aurafix<sup>®</sup>**  
ORTHOPEDIC PRODUCTS





## Woven Epicondylitis Elbow Support

**Ref: 304**

Designed anatomically, manufactured of flexible, sweat preventing woven fabric. Droplet shaped micro-massaging pad applies compression on the muscle groups adherent to the condyles on both sides of the elbow. Applied local compression reduces tendon tension thus the pain.

### Indications:

- 1- Medial or lateral epicondylitis
- 2- Intra articular liquid increase
- 3- Post-op or injury irritations

### Size;

- S** : 22 - 24 cm
- M** : 24 - 26 cm
- L** : 26 - 28 cm
- XL** : 28 - 30 cm
- XXL** : 30 - 32 cm

**aurafix<sup>®</sup>**  
ORTHOPEDIC PRODUCTS



## **Woven Varus Ankle Support**

**Ref: 403**

Manufactured with flexible and sweat-preventing woven fabric. With the elastic bandage, it treats the varus position (stepping with the ankle bent outwards) passively.

### **Size;**

**S** : 19 - 21 cm

**M** : 21 - 23 cm

**L** : 23 - 25 cm

**XL** : 25 - 27 cm

**XXL** : 27 - 29 cm

**aurafix<sup>®</sup>**  
ORTHOPEDIC PRODUCTS



## Air Walking Boot

**Ref: 452**

Intermittent pneumatic compression not only promotes callus formation over fracture site but also reduces local swelling. Shock absorbing sole reduces heel-strike impact force. Front panel for more complete protection and support. Rocker bottom design promotes natural gait. Suitable for wear from contoured boot shell

### Indications:

- 1-** Grade II/III ankle sprain
- 2-** Stable fracture of the fibula, talus, calcaneus and medial malleolus
- 3-** Post ligament, soft tissue and tendon surgery
- 4-** Forefoot or midfoot injury

### Functions:

Aid propulsion  
Substitute ankle motion  
Decreased heel-strike force on calcaneus

**aurafix<sup>®</sup>**  
ORTHOPAEDIC PRODUCTS



## Woven Lumbosacral Corset With Extra Support

**Ref: AO-86**

Manufactured of respirable knitted fabric. With 4 pieces of anatomically designed steel underwire and plastic boost, supports the lumbar region and relieves the stress on the disc thus reduces pain.

### Indications:

- 1-** Lumbago, disc lesions
- 2-** Spondylodesis, spondylolisthesis
- 3-** Post-op preservation and immobilization purposes

### Size;

- S** : 75 - 85 cm
- M** : 85 - 95 cm
- L** : 95 - 105 cm
- XL** : 105 - 115 cm
- XXL** : 115 - 125 cm



**aurafix<sup>®</sup>**  
ORTHOPEDIC PRODUCTS



## **Lumbostad Corset (26cm) With Abdominal Support**

**Ref: AO-61**

Manufactured with flexible, sweat-preventing cotton fabric. In the back, it has 4 pieces of anatomically convenient steel underwire. The front side is composed of two pieces. The lower and thicker piece supports the lower abdomen. The thinner piece on the front side is fastened in a transverse manner. With 26 cm lumbar part, this corset is better choice for shorter patients.

### **Size;**

- S** : 75 - 85 cm
- M** : 85 - 95 cm
- L** : 95 - 105 cm
- XL** : 105 - 115 cm
- XXL** : 115 - 125 cm



**aurafix<sup>®</sup>**  
ORTHOPEDIC PRODUCTS

## Dorsolumbar Corset (26 cm)

**Ref: AO-91**

Manufactured of flexible, respirable woven fabric. Steel underwire on both sides of the vertebral column, extending up to the back and the additional strap with transverse fasteners in the lumbar area provide stabilization. Different from other Dorsolumbar corsets, the axillary strap is made of terry pads which provide a more comfortable usage. The cover and the buckle system ease the wear. The velcro fasteners in the dorsal area provide height adjustment, thus prevent axillary pressure. Terry pads soften the rear part of the corset.

### Indications:

- 1- Vertebral trauma
- 2- Osteoporosis
- 3- Mild kyphosis

### Size;

- S** : 75 - 85 cm
- M** : 85 - 95 cm
- L** : 95 - 105 cm
- XL** : 105 - 115 cm
- XXL** : 115 - 125 cm



**aurafix®**  
ORTHOPEDIC PRODUCTS